

XLMTM_D.2.2.002

Please quote this SOP in your Methods.

Respiratory assessment of a canine model of XLMTM using non-invasive respiratory impedance plethysmography

SOP (ID) Number	XLMTM_D.2.2.002
Version	1
Issued	March 26 th , 2014
Last reviewed	July 10 th , 2017
Author	Melissa A Goddard University of Washington, Seattle WA
Working group members	Martin K Childers, David L Mack University of Washington, Seattle WA
SOP responsible	Melissa A Goddard
Official reviewer	Martin K Childers

TABLE OF CONTENTS

OBJECTIVE.....	3
SCOPE AND APPLICABILITY.....	3
CAUTIONS.....	3
MATERIALS	3
METHODS	4
EVALUATION AND INTERPRETATION OF RESULTS.....	8
REFERENCES	8

1. OBJECTIVE

Anesthetised dogs are fitted with wired bands and connected to a calibrated pneumotachometer at the time of assessment. After a baseline ventilatory pattern is established and measured, the respiratory stimulant doxapram hydrochloride is given. Ventilatory parameters are assessed from measures of airflow and by impedance plethysmography as a function of time in response to doxapram.

2. SCOPE AND APPLICABILITY

Respiratory impedance plethysmography (RIP) and pneumotachography have been used successfully in larger animal models to measure respiration. Similar non-invasive systems have been applied clinically in patients with neuromuscular disorders to assess the respiratory dysfunction that is often present in these diseases.

3. CAUTIONS

The dog should not be shifted or moved during measurement as can cause noise or disruptions in the RIP signal

4. MATERIALS

- EKG pads
- EMKA underjacket and outer jacket of an appropriate size (typically the underjacket is one size smaller; all come in small, medium and large)(EMKA Technologies)
- Alcohol
- Shaver
- Wired band material (EMKA Technologies)
- Scissors
- Abdominal, thoracic and EKG leads (EMKA Technologies)
- Pneumotachometer (EMKA Technologies)
- Calibration unit and t/r box (EMKA Technologies)
- Antennae (EMKA Technologies)
- Small 'screw driver' (included with calibration unit) (EMKA Technologies)
- Laptop, with iox2 software 2.8.0.13 (Windows XP Professional, service pack 2 or 3)
- Telemetry pack (EMKA Technologies)
- 1L super syringe (EMKA Technologies)
- Tubing
- AA Duracell Procell batteries
- Electrode gel (optional)
- Software key (iox) (EMKA Technologies)


XLMTM_D.2.2.002

- Doxapram hydrochloride
- Saline
- Syringe
- Needle
- SD card

5. METHODS

Equipment Set Up

- Ensure that equipment is properly connected and functioning.
 - T/r box connections should be as follows when viewed at the back:



- Lights on the front of the T/r unit should glow green
- Signal strength should be green or high yellow when telemetry box is on
- Connect the antenna
- Turn on the unit
- Connect SD card to T/l box and insert into laptop

XLMTM_D.2.2.002

Load saved config file

- Insert the software key and open the Iox Software
- Load saved configuration file created by installing technician

- Define subject names (Optional)
 - Select **Run** from the main menu to display subject name
 - Input the name or number of the test animal, with a unique modification so there is no duplication in the data directory

- Click **Run!**
- Select **Tuning** → **Calibrate**

Calibrate the inputs

- Calibration values should load automatically as follows for
 - EKG
 - Act
 - RIP (Thor and Ab) bands
 - RIP sum channel

- Uncalibrated inputs will show "--" as their values in the calibration main screen under "Last calibration date/time." If this occurs, values can be set as follows:
 - Select "2 points (manual) calibration"
 - Enter the **Low** and **High values**
 - Enter the unit of measure in **Unit**

	Low	High	Unit
EKG	0	1	mV
Act	0	1	G
	-5V	-1.45V	
RIP	0	1000	m/s
	0	1V	
RIP Sum	0	1000	m/s

- Click **Ok**

XLMTM_D.2.2.002

- Calibrate RF Calc, the pneumotach signal
 - Select **Tuning** → **Calibrate** from the run menu.
 - Select “**integrate calibration**” using the dial
 - Enter the **Low value** (0) and the **High value** (1000)
 - Enter the unit of measure in **Unit** (ml/sec or l/sec)
 - Set the calibrating unit/amplifier to **DC** using the switch
 - Set the gain to **1**
 - Adjust the coarse adjustment (lowest screw) using the provided ‘screw driver’ until the signal is close to **0V** (no flow)
 - Adjust the gain in a stepwise fashion, dialing back until the signal returns to 0 using the fine adjustment dial
 - Set the **integration period** to the number of seconds sufficient to inject the volume. 2s or 3s is ideal for accurate calibration, although longer periods can be used initially.
 - Attach the 1L syringe to the pneumotachometer using the tubing. When connected correctly, pushing air in should cause a negative deflection to be measured
 - Flip the metal transducer if this is not observed and reconnect the pneumotach
 - Inject the volume and immediately click the **High** button
 - Ensure that the signal is within range and that injection time does not exceed the integration period. (Area under the curve should be between the 3s.)
 - If not, dial up the gain until the scale fits.
 - Switch the calibration unit back to **AC** and deselect “remove AC offset” on the program
 - Click **Ok**
- Save this new configuration file under a new name (which includes the date of creation)

Preparation of Dogs and Jacket Fitting

- Shave the anaesthetised dog as cleanly as possible in the appropriate areas :
 - under the left ‘armpit’ of the forelegs
 - under the right ‘armpit’ of the forelegs
 - On the left side of the groin
- Size the dog for the thoracic and abdominal bands in such a way that they fit snugly but still allow two fingers under the band
 - Cut using the scissors

XLMTM_D.2.2.002

- Wipe the area with alcohol to remove any stray hairs and fix the EKG pads firmly to the dog. Electrode gel can also be used to improve adhesion and conductivity.
- Connect the EKG leads as follows:
 - Green** at the **groin**
 - White** on the **right**
 - Black** on the left
- Feed the band material through the loops of the inner jacket and place the inner jacket on the dog, so that the bands face outwards. The metal part of the bands should also be visible (ie, facing outwards)
- Check to ensure that the thoracic and abdominal leads are in an appropriate 'open' position, feed the end of the band material into each and then close. Ends should lock firmly, but should not need to be forced as this may cause damage
 - The blue measures the thoracic band
 - The yellow measures the abdominal
 - Dog should appear as follows:



- Place the outer jacket on the dog (Optional)
- Connect thoracic, abdominal and EKG leads to the telemetry pack. Connecting the EKG lead turns the unit on.
- Ensure that a signal is being obtained on the computer
- Place the telemetry unit in the pouch of the outer jacket before continuing with any planned experiments (Optional)
- A catheter should also be placed in the fore or hindlimb for the injection of doxapram hydrochloride.

XLMTM_D.2.2.002

Run the experiment and start streaming and storage

- To start **streaming**, select **Start/Stop → Storage/Streaming**,
 - Click **start streaming**
 - Check the **Save beats** box
 - Click **Ok**
- Or click the pre-programmed (**F1**) key
- To stop streaming press **F2**

- To start **storage**, select **Start/Stop → Storage/Streaming**,
 - Click **start storage**
 - Select the site or press **Select All sites**
 - Click **Ok**
- Or click the pre-programmed (**F3**) key
- To stop storage press **F4**

Doxapram hydrochloride administration

- Create a label in Iox by pressing the **F12** button and entering the amount of doxapram.
- Administer a bolus injection of 1.0mg/kg of doxapram hydrochloride, followed rapidly by a saline flush
 - Press **Enter** at the same time the dose is administered to label
- Collect a minimum of 20 clear breaths once respiration has reached its peak

6. EVALUATION AND INTERPRETATION OF RESULTS

Iox and ecgGAUTO (version 2.8.2.6) (EMKA Technologies) calculates various respiratory parameters from the measured curves such as respiratory rate, tidal volume, peak inspiratory and expiratory flows, and thoracic and abdominal volume. The software can also be used to generate reports of these values, including flow-volume or phase curves.

7. REFERENCES

Iox2 user manual. *EMKA technologies*. V2.4. rev 2. 2010

Goddard MA, Mitchell EL, Smith BK, Childers MK. Establishing clinical end points of respiratory function in large animals for clinical translation. *Phys Med Rehabil Clin N Am*. 2012;23(1):75-94, xi.